



**TOGETHER
TOWARDS
HEALTH**

**ASSESSING AND PREPARING FOR FUTURE
REQUIREMENTS AND STANDARDS OF
PHARMACEUTICAL SOLUTION BAGS**

ABRASP Conference December 7, 2021



Rely on it.



LET'S SHAPE THE FUTURE OF HEALTHCARE



Peter Robben
Global Product Manager
RENOLIT Healthcare



Diego Morais
General Manager
RENOLIT Chile SpA



Rely on it.



WE ARE PART OF THE RENOLIT GROUP

The **RENOLIT Group**
stands for partnership-based
service, an international focus and
innovative, high-quality, reliable products.



Rely on it.



OUR NEW STRATEGY AND ITS FOUR CORNER STONES

COOPERATION



INNOVATION



SUSTAINABILITY



WORLDWIDE



Rely on it.



LOOK NO FURTHER – WE ARE YOUR GLOBAL SPARRING PARTNER



1 RENOLIT SE – HEADQUARTERS
and RENOLIT SE (Thansau)

2 RENOLIT Nederland B.V. and
RENOLIT Amsterdam Innovation Center

3 American RENOLIT Corporation – LA and
RENOLIT California Innovation Center

4 RENOLIT Chile SpA

5 RENOLIT Healthcare Beijing and
RENOLIT Beijing Innovation Center



Rely on it.



OUR MOST COMPREHENSIVE HEALTHCARE PORTFOLIO EVER

APPLICATION FIELDS



Blood & blood components



Biotechnology



Nutrition



IV & pharmacological applications



Dialysis




Further exciting areas are in development

COMPLEMENTARY



Ports & caps




Barrier films & secondary packaging


SERVICES



Custom-made solutions



Tailor-made films



EMT equipment lease



Rely on it.

FLEXIBLE BAGS FOR PHARMACEUTICAL SOLUTIONS

OUTLINE

- » **Overview** of the polymeric materials typically used.
- » **Guidance** for selection of polymeric components for use in bags:
 - » Film
 - » Tubing
 - » Ports
- » **Requirements** for suppliers of polymeric components used in the construction of bags for pharmaceutical solutions.





TOGETHER
TOWARDS
HEALTH



MATERIAL SELECTION CRITERIA

MATERIAL SELECTION CRITERIA

BAGS FOR PHARMACEUTICAL SOLUTIONS

- » Typically comprised of film, ports and / or tubing.
- » Subject to multiple steps:
 - » Manufacturing of components (film, tubing, ports, label, ...)
 - » Assembly (bag sealing)
 - » Filling
 - » Closing
 - » Overwrapping
 - » Sterilization
 - » Transport & storage
- » Subject to multiple requirements:
 - » Shelf life
 - » Drug compatibility
 - » Cleanliness
 - » Environmental
 - » In-use performance
 - » Validation
 - » Change Control
 - » ...



MATERIAL SELECTION CRITERIA

MANY OPTIONS FOR SELECTION OF COMPONENTS FOR FLEXIBLE BAGS

- » Broad range of material options available
- » Different performance and properties
- » Shape/form
- » Multilayer structure of film and tubing
- » Port design

COMMON POLYMERS USED IN FLEXIBLE BAGS FOR PHARMACEUTICAL SOLUTIONS

- | | |
|-------------|------------------------|
| » PP | » PET |
| » PE | » Tie layer copolymers |
| » EVA | » ... |
| » PVC | |
| » PVDF, FEP | |
| » EVOH | |
| » PA | |
| » SEBS | |





MATERIAL SELECTION CRITERIA

KEY MATERIAL-RELATED PROPERTIES

- » Sterilization resistance: steam, radiation or ETO
- » Operational temperature
- » API compatibility
- » Gas permeability (WVTR, OTR, CO₂TR)
- » Bag manufacturing technology



Rely on it.

MATERIAL SELECTION CRITERIA

	PVC ¹	PP ²	EVA ³	PE
Operational temperature	-30° – 90°C	-20° – 100°C	-196° – 80°C	-80° – 80°C
Sealing technology	HF	Heat FFS	HF	Heat
Sterilization of filled bags	Steam	Steam	–	–
Sterilization of empty bags	–	–	EtO, Radiation	Radiation
Fluid bags	Terminal steam sterilization	Terminal steam sterilization	Aseptic filling	Aseptic filling
Typical pharmaceutical use	IV, Dialysis	IV, Dialysis, Nutrition	Nutrition	Biotech SUS
API compatibility	-- to ++	++	0 to ++	++
M/C bags peelable seal	no	yes	no	no
WVTR	high	medium	high	medium
OTR	high	high	high	high

¹ Plasticizer selection influences PVC properties: gas permeability, hardness / flexibility, appearance (optical, odor), fluid compatibility / E&L


² PP films are typically multi-layer structures and including other polymers

³ EVA properties vary with copolymer content

MATERIAL SELECTION CRITERIA

FILM


- » Multilayer structure
 - » Combination of different polymers
 - » Tailored properties/performance
 - » Peelable seals for multi-chamber bags
 - » O₂/CO₂/UV barrier
- » Double wound, single wound or tubular,
 - » Different bag manufacturing equipment
 - » Different film production processes
- » Typical film thickness for fluid bags
 - » PP 0.18 – 0.25 mm
 - » EVA 0.30 – 0.40 mm
 - » PVC 0.25 – 0.40 mm



MATERIAL SELECTION CRITERIA

PORT SYSTEM

- » Tubing with ports
 - » Port assembled in tubing
 - » Adhesive bonding PC to PVC
 - » Heat bonding during steam sterilization
 - » PC to PVC
 - » PC to PP multilayer tube
 - » PP to PP mono/multilayer tube
 - » Option for variety in ports
- » Bag port system
 - » Port sealed to the film
 - » Port closed (sealed) after filling
 - » Less variety in ports; bag making machine modification
 - » In general lower cost





TOGETHER
TOWARDS
HEALTH



SUPPLIER SELECTION CRITERIA



SUPPLIER SELECTION CRITERIA: QUALITY MANAGEMENT

THERE IS A LOT MORE THAN MATERIAL PERFORMANCE AND PRICING...

- » Pharma-oriented Quality Management System
 - » ISO 13485, ISO 15378 certification
 - » Periodic customer audits
 - » Validation process
 - » Management of changes
 - » Customer notifications
 - » Management of non-conformances
 - » Traceability
 - » CoC / CoAs
 - »
- » Cleanliness
 - » Component surface exposed to environment during production
 - » ISO class 7 cleanroom production minimum preferred.



SUPPLIER SELECTION CRITERIA: REGULATORY SUPPORT

THERE IS A LOT MORE THAN MATERIAL PERFORMANCE AND PRICING...

- » Regulatory support
 - » DMF
 - » Statements on (absence of) substances of concern
 - » Pharmacopoeia compliance
 - »
- » Technical support
 - » Processing, performance troubleshooting
 - » Relevant biocompatibility reports
 - » ISO 10993
 - » Ph.Eur., JP, USP, ...
 - » Physical/chemical properties
 - » E&L
 - » USP 661.1
 - » BPOG (Single-use bioprocessing)






SUPPLIER SELECTION CRITERIA: SECURITY OF SUPPLY

THERE IS A LOT MORE THAN MATERIAL PERFORMANCE AND PRICING...

- » Security of supply
 - » IT HAPPENS
 - » Backup options
 - » Raw materials
 - » Component(s)
- » Multiple production sites
 - » Typically covers extra set of regional raw materials
 - » Supports global manufacturing strategy



Rely on it.

SUPPLIER SELECTION CRITERIA: CHANGE MANAGEMENT

THERE IS A LOT MORE THAN MATERIAL PERFORMANCE AND PRICING...


- » Consistent, reliable supply is ideal. BUT change happens.
- » Supplier Change Management is fundamental
 - » Change control and notification process in place
 - » Timeliness of notification
 - » Last buy options
 - » Alternative products ready to go
- » Supplier understanding of the pharmaceutical industry is essential
- » Transparency/open communication always helps
- » BPSA/BPOG Change Notification guide for Single-Use Bioprocessing Systems



MATERIAL AND SUPPLIER SELECTION CRITERIA

SUMMARY

- » Pharmaceutical Solution Bags are not simple
 - » Many material options available
 - » It's more than a piece of plastic
 - » Changes can be long and complex
- ▶ Conduct a broad assessment of material and supplier requirements.
 - ▶ Engage with suppliers early in the product development process
 - ▶ Suppliers with a strong pharmaceutical experience and pharma industry business orientation are ideal.
 - ▶ Design with supply security embedded from the beginning.





**TOGETHER
TOWARDS
HEALTH**



**TOGETHER WE TAKE
HEALTHCARE TO THE NEXT LEVEL**

THANK YOU FOR YOUR ATTENTION



Rely on it.